

Lung-RADS v2022 Clinical Guidance Module

Executive Summary

The Lung-RADS v2022 Clinical Guidance module is a structured decision support application embedded within Nuance PowerScribe. It provides radiologists with real-time, guideline-concordant Lung-RADS category assignment and management recommendations for lung cancer screening CT examinations. The module encodes the complete decision logic from the American College of Radiology's Lung-RADS v2022 guideline.

The module operates within the PowerScribe Clinical Guidance framework. It supports input via dropdown selectors and voice entry, with each selection mapped to a discrete lexicon element defined by the guideline. The module dynamically resolves the appropriate Lung-RADS category (0–4X) based on encoded rule pathways and returns the corresponding management recommendation. Each rule is directly derived from published ACR consensus criteria.

The module inserts standardized category and recommendation text into the report's Impression section, formatted for readability and accompanied by a citation to the source publication.

This document defines the design, logic structure, input specification, output schema, and interaction model for the Lung-RADS v2022 module. It is intended for implementation stakeholders, including clinical informaticists, software engineers, QA specialists, and reviewers.

Clinical and Technical Scope

Clinical Scope

The module encodes the complete classification and management framework defined in the ACR Lung-RADS v2022 guideline. It applies exclusively to low-dose chest CT performed for lung cancer screening in high-risk patients; it is not intended for incidental pulmonary nodule reporting.

Lesion types included in scope:

- Solid nodules (baseline, new, growing, or stable)
- Non-solid (ground-glass) nodules
- Part-solid nodules with defined solid components
- Atypical pulmonary cysts (thin-walled, thick-walled, multilocular, or progressive)
- Juxtapleural/perifissural nodules with triangular/oval/lentiform morphology
- Airway nodules (subsegmental, segmental, or proximal)
- Special handling for stable/decreasing nodules, proven benign lesions, and slow-growing nodules
- Application of the “S” modifier for clinically significant incidental findings
- Application of “4X” for nodules with additional suspicious features

The module also encodes management recommendations ranging from annual screening to short-interval follow-up CT to diagnostic work-up (diagnostic CT, PET/CT, biopsy, or surgical referral).

The logic conforms as exactly as possible to the published guideline. No local modifications or institution-specific heuristics are applied.

Technical Scope

The module is implemented as a Clinical Guidance (CG) module within Nuance PowerScribe. It is accessible during report creation and can be invoked by voice or macro trigger. Input is collected through a structured interface using dropdown menus, each linked to a defined lexicon value. Voice navigation is supported for all selectable fields.

The output is generated as standardized text and inserted into the Impression section, including an embedded citation to the Lung-RADS v2022 publication. The impression content is dynamically assembled based on selected features and resulting category, with grammatical variation to support narrative readability.

Out of scope for this module:

- Extraction of nodule characteristics from free-text.
- Application to non-screening CT studies.
- Customization of management recommendations beyond guideline-defined pathways.

Guideline Background and Source Authority

The Lung-RADS Clinical Guidance module is derived from the American College of Radiology's Lung Imaging Reporting and Data System (Lung-RADS®) v2022. This guideline defines a standardized lexicon, decision framework, and management schema for classification of findings on low-dose CT performed for lung cancer screening.

Purpose of the Source Guideline

Lung-RADS was developed to:

- Standardize CT terminology for describing pulmonary nodules in the screening setting.
- Assign numeric categories (0–4X) based on nodule size, growth, morphology, and context.
- Link each category to a management recommendation that balances early detection with avoidance of unnecessary testing.

The v2022 update incorporates:

- Clarified definitions for airway nodules.
- Formal criteria for atypical pulmonary cysts, including wall thickness and multilocularity.
- Explicit guidance for slow-growing nodules.
- Clarifications on growth definition (>1.5 mm in 12 months) and size threshold crossing.
- Updated recommendations for use of the “S” modifier for significant incidental findings.

Implementation Scope within the Module

This module encodes the guideline in full, without omission or modification. The decision logic reflects all assessment categories:

- 0 – Incomplete
- 1 – Negative
- 2 – Benign

- 3 – Probably Benign
- 4A – Suspicious
- 4B – Very Suspicious
- 4X – Suspicious with Additional Features

The module supports classification of solid, non-solid, part-solid, airway, juxta-pleural, and cystic lesions, with downgrade rules for stable/decreasing nodules and reversion for proven benign findings.

Citation

ACR Lung-RADS v2022: Assessment categories and management recommendations. Christensen J, Prosper AE, Wu CC, et al. J Am Coll Radiol. 2024;21(3):473–488. doi:10.1016/j.jacr.2023.09.009

Available online: [https://www.jacr.org/article/S1546-1440\(23\)00761-5/fulltext](https://www.jacr.org/article/S1546-1440(23)00761-5/fulltext)

Data Inputs

The Lung-RADS module accepts a fixed set of structured inputs corresponding to guideline descriptors. All inputs are collected via dropdown menus or voice-selectable fields.

Required Inputs

- Lesion type (solid, non-solid, part-solid, cystic, juxta-pleural, airway)
- Nodule status (baseline, new, growing, stable, decreasing)
- Mean nodule diameter (mm)

Optional Inputs

- Lobe location
- “4X” upgrade flag
- “S” modifier for incidental findings
- Validation ensures that all required fields are present before category assignment.

The screenshot shows a web-based input form for Lung-RADS v2022. It features a light gray background with a dark blue border. On the left, there are labels for required and optional inputs. On the right, there are dropdown menus and a list of options. The 'Dominant Lesion Type' dropdown is set to 'Part Solid'. The 'Nodule Status' dropdown is open, showing a list of options: 'Initial/Baseline', 'New', 'Stable Overall, Stable Solid Component', 'Stable Overall, Growing Solid Component', 'Growing Overall, Stable Solid Component', 'Growing Overall, Growing Solid Component', 'Slow-Growing (<1.5mm/yr over mult. exam)', and 'Decreasing in size'. The 'Mean diameter (mm)' and 'Mean solid diameter (mm)' fields are empty. The 'Optional fields' section is also empty.

✓ Dominant Lesion Type:	Part Solid
Nodule Status:	Initial/Baseline
Mean diameter (mm):	New
Mean solid diameter (mm):	Stable Overall, Stable Solid Component
	Stable Overall, Growing Solid Component
	Growing Overall, Stable Solid Component
	Growing Overall, Growing Solid Component
	Slow-Growing (<1.5mm/yr over mult. exam)
Optional fields:	Decreasing in size

Figure 1. The module comprises required and optional inputs to drive generation of impression text.

Decision Logic and Category Assignment

The module encodes the full categorical logic of Lung-RADS v2022. Logic execution is deterministic, and each terminal pathway resolves to a single category with a linked management recommendation.

Decision Pathways:

- The most suspicious nodule determines category.
- Growth: >1.5 mm increase in mean diameter within 12 months; slow growth over years → 4B.
- Threshold crossing: size change sufficient to enter a higher-risk category prompts upgrade.
- Stable/decreasing nodules: downgraded (4A→3, 3→2).
- Airway nodules: subsegmental = 2; segmental/proximal = 4A, persistent = 4B.
- Atypical cysts: thin-walled = benign; thick-walled/multilocular or enlarging = 4A/4B.
- Infection/inflammation: findings suspicious for acute process may be assigned Category 0 with 1–3 month follow-up.

Output Generation

The Lung-RADS module inserts standardized content into the Impression section of the radiology report.

Output components include:

- Category statement (e.g., “Lung-RADS Category 3: Probably benign.”)
- Lesion description (nodule type, size, stability)
- Management recommendation (e.g., “Follow-up LDCT in 6 months.”)
- Citation

If required inputs are missing, no output is inserted and the user is prompted to review.

The screenshot shows the 'Guidance: Lung RADS 2022' window. On the left, there are input fields for 'Quick Select (Lung-RADS Chart)', 'Dominant Lesion Type' (set to 'Part Solid'), 'Nodule Status' (set to 'Initial/Baseline'), 'Mean diameter (mm)' (set to '7'), and 'Mean solid diameter (mm)' (set to '4'). Below these are 'Optional fields' including 'Dominant nodule lobe' (set to 'RLL'), 'Upgrade to Category 4B?' (set to 'No'), and 'Significant findings?'. On the right, the 'Include' dropdown is set to 'Findings, Impression, Recommendation'. The 'Impression' section displays: 'Lung-RADS Category: 3, probably benign based on imaging features or behavior. The dominant lesion is a part-solid nodule measuring 7 mm with a solid component measuring 4 mm located in the right lower lobe.' The 'Recommendation' section displays: 'Follow-up Lung Cancer Screening Chest CT in 6 months.'

Figure 2. Example impression and recommendation text.

User Interface and Interaction Behavior

- Invoked by voice command, sidebar selection, or macro.
- Structured input panel guides completion of required and optional fields.
- Logic executes automatically once required fields are present.
- Quick Select mode allows direct category assignment.
- An optional interactive flowchart mirrors the Lung-RADS decision structure, allowing visual navigation and direct category selection.

Dominant Lesion Type:

Optional fields:

Significant findings?

Juxtapleural (<10mm, solid, smooth)

Juxtapleural (<10mm, solid, smooth)

Cat 4B Proven Benign

Cat 3/4A Resolved on Follow-up

Atypical Pulmonary Cyst

Airway Nodule

Nodule w/Benign Features

Infection/Incomplete/Priors

No Nodules Present

Nuance PowerScribe One

File Edit View Insert Format Tools Speech Help

Close Draft Correct Reject Prelim Sign Discard

AutoText New...

Report - UNKNOWN, - anvrvjkl

Ent

Guidance: Lung RADS 2022

Quick Select (Lung-RADS Chart):

:OR:

Complete the fields below:

✓ Dominant Lesion Type: Solid

Nodule Status:

Mean diameter (mm):

Optional fields:

Dominant nodule lobe:

Significant findings?

Clinical Guidance

Lung-RADS	Category Descriptor	Findings
0	Incomplete Estimated Population Prevalence: ~1%	Prior chest CT examination being located for comparison (see note 9) Part or all of lungs cannot be evaluated Findings suggestive of an inflammatory or infectious process (see note 10)
1	Negative Estimated Population Prevalence: 39%	No lung nodules OR Nodule with benign features: - Complete, central, popcorn, or concentric ring calcifications OR - Fat-containing
2	Benign Based on imaging features or indolent behavior Estimated Population Prevalence: 45%	Juxtapleural nodule: - < 10 mm (< 524 mm ³) mean diameter at baseline or new AND - Solid; smooth margins; and oval, lentiform, or triangular shape Solid nodule: - < 6 mm (< 113 mm ³) at baseline OR - New < 4 mm (< 34 mm ³) Part-solid nodule: - < 6 mm total mean diameter (< 113 mm ³) at baseline Non-solid nodule (GGN): - < 30 mm (< 14,137 mm ³) at baseline, new, or growing OR - ≥ 30 mm (≥ 14,137 mm ³) stable or slow-growing (see note 7) Airway nodule, subsegmental at baseline, new, or stable (see note 11) Category 3 nodule that is stable or decreased in size at 6-month follow-up CT, OR Category 3 or 4A nodules that resolve on follow-up, OR Category 4B findings proven to be benign in etiology following appropriate diagnostic workup
3	Probably Benign Based on imaging features or behavior Estimated Population Prevalence: 9%	Solid nodule: - ≥ 6 to < 8 mm (≥ 113 to < 268 mm ³) at baseline OR - New 4 mm to < 6 mm (34 to < 113 mm ³) Part-solid nodule: - ≥ 6 mm total mean diameter (≥ 113 mm ³) with solid component < 6 mm (< 113 mm ³) at baseline OR - New < 6 mm total mean diameter (< 113 mm ³) Non-solid nodule (GGN): - ≥ 30 mm (≥ 14,137 mm ³) at baseline or new Atypical pulmonary cyst: (see note 12) - Growing cystic component (mean diameter) of a thick-walled cyst Category 4A nodule that is stable or decreased in size at 3-month follow-up CT (excluding any airway nodules) Solid nodule:

Figures 3 and 4. The user may either input nodule information by field or select a clickable Lung-RADS chart to prepopulate fields.

